

## **GIT 4 STOMACH WITH DOUBLE LIGATION**

## **SALIENT FEATURES:**

- Two ligatures are applied at the **cardiac end** of the esophagus and two ligatures below the **pyloric end** of the stomach.
- The stomach is removed by cutting between the double ligatures at both ends, and is **opened along the greater curvature**.
- The **mucous membrane** is examined for the presence of any stain, congestion, hemorrhage, desquamation, ulceration, sloughing or perforation.
- The content of the stomach is noted in respect to quantity, nature of material/food, state of digestion, color and smell.
- STOMACH: length 25 to 30 cm. Weight: 300 g; capacity 1100 to 1200 ml